


MIRADO™

Hydrophobic Foldable Intraocular Lens

TM

Hydrophobic Lens :

- A hydrophobic lens without compromise for mini-incisions
- Excellent material & design which allows smooth implantation
- Easy implantation – 2.2-2.4 mm incision
- Compatible with Fully-Preloaded Implantation System at 2.2-2.4 mm
- It is a safe and reliable IOL resulting in efficient minimal invasive surgery
- Controlled unfolding of haptics and excellent centration
- The square edge of the lens helps to minimize PCO



Optic Design - Aspheric/Bi-Convex,
Edge Design - 360° Continuous
Square Edge, Haptic Angle - 0°/5°, °
Diopter Range : +6D to +28.5D(15D
to 28D in 0.5D increments)
*Available In Clear & Natural Yellow
**Other Diopters Available On
Request


Also Available in Preloaded System

MODEL NO.-HLL260125CRCA

C-Loop

KEY SPECIFICATIONS

Optic Diameter	6.00 / 6.50 mm
Overall Diameter	12.00 / 12.50 / 13.00 mm
Haptic Design	C-Loop




MODEL NO.-HCL160130CRCA

L-Loop

KEY SPECIFICATIONS

Optic Diameter	6.00 mm
Overall Diameter	12.00 / 12.50 / 13.00 mm
Haptic Design	L-Loop



OPTIMIZING YOUR SURGICAL RESULT

MIRADO BIOMETRY VALUES

A Constant (Ultrasound technique)	117.8
A Constant (Optical technique) SRK/T	118.8
Holladay 1 –sf	1.66
Hoffer Q –pACD	5.45
Haigis	a0=1.26 a1=0.41 a2=0.10

