


CANNULAS



- ANESTHESIA
- ANTERIOR CHAMBER IRRIGATION
- CAPSULE POLISHING
- CORTEX REMOVAL
- CAPSULOTOMY
- HYDRODISSECTION
- IRRIGATION AND AIR INJECTION
- NUCLEUS REMOVAL
- IRRIGATION AND ASPIRATION
- ANTERIOR CHAMBER MAINTENANCE


ANESTHESIA

RETROBULBAR, Atkinson Tip





Code	Description	Contents
111638	23G x 1½ in. (0.6 x 38mm)	10/Box
111637	25G x 1½ in. (0.5 x 38mm)	10/Box

RETROBULBAR, Sharp Tip





Code	Description	Contents
111278	25G x 1 in. (0.5 x 25mm)	10/Box
111282	30G x ½ in. (0.3 x 12.5mm)	10/Box

PERIBULBAR, Atkinson Tip



Code	Description	Contents
115144	23G x 7/8 in. (0.6 x 22mm)	10/Box
115027	25G x 7/8 in. (0.5 x 22mm)	10/Box
111642	25G x 1¼ in. (0.5 x 32mm)	10/Box


SUB-TENON



Code	Description	Contents
125176	19G x 7/8 in. (1.1 x 22mm), open tip	10/Box



ANTERIOR CHAMBER IRRIGATION

ANGLED TIP, 4mm (Rycroft)





Code	Description	Contents
121281	20G (0.9mm)	10/Box
121639	25G (0.5mm)	10/Box
101280	27G (0.4mm)	10/Box
121273	30G (0.3mm)	10/Box

ANGLED TIP, 6mm (Rycroft)





Code	Description	Contents
121287	20G (0.9mm)	10/Box
121286	27G (0.4mm)	10/Box
121276	30G (0.3mm)	10/Box

ANGLED TIP, 7mm (Rycroft)




Code	Description	Contents
121637	25G (0.5mm)	10/Box
121290	27G (0.4mm)	10/Box
121277	30G (0.3mm)	10/Box

ANGLED TIP, 8mm (LASIK)



Code	Description	Contents
121293	27G (0.4mm)	10/Box



ANGLED TIP, 10mm (LASIK)



Code	Description	Contents
121280	25G (0.5mm)	10/Box
121274	27G (0.4mm)	10/Box



CAPSULE POLISHING

SILICONE TIP, Angled 7mm




Code	Description	Contents
175041	23G (0.6mm)	10/Box
175005	27G (0.4mm)	10/Box

POSTERIOR, Angled 9mm, Top Port





Code	Description	Contents
171604	21G (0.8mm)	10/Box

POSTERIOR, Angled 7mm, Open Tip



Code	Description	Contents
171614	21G (0.8mm)	10/Box


OLIVE TIP



Code	Description	Contents
171589	23G (0.6mm)	10/Box


CAPSULOTOMY

FORMED





Code	Description	Contents
131599	23G x $\frac{5}{8}$ in. (0.6 x 16mm)	10/Box
131610	25G x $\frac{5}{8}$ in. (0.5 x 16mm)	10/Box
131618	27G x $\frac{5}{8}$ in. (0.4 x 16mm)	10/Box

ANGLED 45°




Code	Description	Contents
135049	25G x $\frac{5}{8}$ in. (0.5 x 16mm)	10/Box
135093	27G x $\frac{5}{8}$ in. (0.4 x 16mm)	10/Box

REVERSE BEND




Code	Description	Contents
135115	25G x $\frac{5}{8}$ in. (0.5 x 16mm)	10/Box

SHORT





Code	Description	Contents
135004	25G x $\frac{5}{8}$ in. (0.5 x 16mm)	10/Box
135008	27G x $\frac{5}{8}$ in. (0.4 x 16mm)	10/Box

SIDE-CUT



Code	Description	Contents
131602	25G x $\frac{5}{8}$ in. (0.5 x 16mm)	10/Box
131609	27G x $\frac{5}{8}$ in. (0.4 x 16mm)	10/Box



PUSH-PULL, Formed



Code	Description	Contents
135042	27G x $\frac{5}{8}$ in. (0.4 x 16mm)	10/Box



HYDRODISSECTION

FLAT TIP, Angled 7mm





Code	Description	Contents
125024	23G (0.6mm)	10/Box
125036	25G (0.5mm)	10/Box
125162	27G (0.4mm)	10/Box

FLAT TIP, Angled 8mm





Code	Description	Contents
125037	25G (0.5mm)	10/Box
125158	27G (0.4mm)	10/Box
121592	30G (0.3mm)	10/Box

HORIZONTAL FLAT TIP, Curved




Code	Description	Contents
125155	27G (0.4mm)	10/Box

VERTICAL FLAT TIP, Curved



Code	Description	Contents
125099	27G (0.4mm)	10/Box



J-SHAPED, Vertical Flat Tip



Code	Description	Contents
125029	27G (0.4mm)	10/Box



HYDRODISSECTION

ANGLED STANDARD TIP (Chang)





Code	Description	Contents
121198	27G (0.4mm), 1.5mm tip	10/Box

CORTICAL CLEAVING, Flat Tip, Angled 11mm





Code	Description	Contents
105158	27G (0.4mm)	10/Box

CORTICAL CLEAVING, Round Tip, Angled 11mm





Code	Description	Contents
125157	27G (0.4mm)	10/Box

CORTICAL CLEAVING, Round Tip, Angled 12mm



Code	Description	Contents
125197	30G (0.3mm)	10/Box



TAPERED, Angled 1mm



Code	Description	Contents
125159	25 to 30G (0.5 to 0.3mm)	10/Box



IRRIGATION AND AIR INJECTION

VISCOELASTIC, Angled





Code	Description	Contents
125009	23G x $\frac{7}{8}$ in. (0.6 x 22mm)	10/Box
125026	27G x $\frac{7}{8}$ in. (0.4 x 22mm)	10/Box

VISCOELASTIC, Formed




Code	Description	Contents
125007	27G x $\frac{7}{8}$ in. (0.4 x 22mm)	10/Box

WASHOUT, Angled 8mm (Bishop-Harmon)



Code	Description	Contents
165040	19G x $\frac{5}{16}$ in. (1.1 x 8mm)	10/Box

OLIVE TIP, Curved



Code	Description	Contents
165147	23G x $\frac{7}{8}$ in. (0.6 x 22mm)	10/Box

NUCLEUS REMOVAL



VECTIS, 1 Port



Code	Description	Contents
141619	24G (0.55mm)	10/Box

IRRIGATION AND ASPIRATION



SIMCOE, Standard



Code	Description	Contents
161571	23G (0.6mm), disposable	10/Box
161573	23G (0.6mm), reusable	1/Box

ANTERIOR CHAMBER MAINTENANCE

RIBBED OPEN TIP



Code	Description	Contents
125061	20G (0.9mm)	3/Box
125062	23G (0.6mm)	3/Box